

Anti-Cd105 Rabbit Polyclonal Antibody For Wb, Elisa - Q63961

Item Number: CM0022507

Introduction

This rabbit polyclonal antibody specifically detects mouse CD105 (Endoglin, S-endoglin). Validated for Western blot (WB) and ELISA, it targets a coreceptor in TGF-beta/BMP signaling and a key regulator of angiogenesis, making it suitable for vascular biology, cancer, and developmental research.

[Learn More](#)

Parameter	Value
Product Name	CD105 Rabbit pAb
Remarks/Alias	Endo; CD105; S-endoglin
Species	Mouse
Gene ID	13805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-581 of mouse CD105 (NP_031958.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q63961
Protein Weight	70kDa
Shipping	Ice bag