

Anti-Cd105/Endoglin Monoclonal Antibody For Wb, Ihc-P, Elisa - P17813

Item Number: CM0002126



Introduction

Recombinant rabbit monoclonal anti-CD105/Endoglin antibody validated for human endoglin detection in WB, IHC-P, and ELISA. Ideal for research on angiogenesis, TGF-beta signaling, and vascular biology.

[Learn More](#)

Parameter	Value
Product Name	CD105/Endoglin Rabbit mAb
Remarks/Alias	END; HHT1; ORW1; CD105/Endoglin
Species	Human
Gene ID (Human)	2022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CD105/Endoglin (P17813).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P17813
Protein Weight	71kDa
Shipping	Ice bag