

Anti-Cd105/Endoglin Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q63961

Item Number: CM0006319

Introduction

High-specificity rabbit monoclonal antibody against mouse CD105/Endoglin, validated for WB, FC, and ELISA. Ideal for angiogenesis research, vascular biology studies, and TGF-beta signaling pathway analysis.

[Learn More](#)

Parameter	Value
Product Name	CD105/Endoglin Rabbit mAb
Remarks/Alias	Endo; CD105; S-endoglin
Species	Mouse
GenelD (Human)	2022
GenelD	13805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-581 of mouse CD105 (NP_031958.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse
SWISS	Q63961
Protein Weight	70kDa
Shipping	Ice bag