

Anti-Cd105/Endoglin Rabbit Polyclonal Antibody For Wb And Elisa - P17813

Item Number: CM0029082



Introduction

Rabbit polyclonal antibody against human CD105/Endoglin, validated for WB and ELISA, cross-reacts with mouse. Suitable for angiogenesis research and endothelial cell studies.

[Learn More](#)

Parameter	Value
Product Name	CD105/Endoglin Rabbit pAb
Remarks/Alias	END; HHT1; ORW1; CD105/Endoglin
Species	Human
Gene ID (Human)	2022
Gene ID	2022
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-337 of human CD105/Endoglin (NP_001108225.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17813
Protein Weight	71kDa
Shipping	Ice bag