

Anti-Cd107B (Lamp2) Rabbit Polyclonal Antibody For Wb, Elisa - P17047

Item Number: CM0023783



Introduction

Rabbit polyclonal antibody targeting mouse CD107b (LAMP2) at ~46 kDa. Validated for WB and ELISA. Suitable for research on lysosomal biogenesis, autophagy, and chaperone-mediated autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	CD107b Rabbit pAb
Remarks/Alias	Mac3, LGP-B, CD107b, Lamp-2, Lamp II, Lamp-2a, Lamp-2b, Lamp-2c
Species	Mouse
Gene ID (Human)	3920
Gene ID	16784
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-189 of mouse CD107b (NP_001017959.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P17047
Protein Weight	46kDa
Shipping	Ice bag