


Anti-Cd109 Monoclonal Antibody For Wb, Elisa - Q6Yhk3

Item Number: CM0007329



Introduction

This rabbit monoclonal antibody targets human CD109, a GPI-anchored cell membrane protein, for Western blot and ELISA. Suitable for investigating TGF-beta signaling modulation in keratinocytes, endothelial cells, and platelets.

[Learn More](#)

Parameter	Value
Product Name	CD109 Rabbit mAb
Remarks/Alias	CD109; CPAMD7; p180; r150; CD109 antigen
Species	Human
GeneID (Human)	1024
GeneID	135228
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 223-462 of human CD109. (NP_598000.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6YHK3
Protein Weight	74kDa/152kDa/159kDa/161kDa
Shipping	Ice bag