

Anti-Cd10/Mme Monoclonal Antibody For Wb, Ihc-P, Elisa - P08473

Item Number: CM0004558



Introduction

Rabbit monoclonal antibody against human CD10/MME (Neprilysin). Validated for Western blot, IHC-P, and ELISA. Suitable for studying neuropeptide degradation, cardiovascular peptide inactivation, and tumor antigen research.

[Learn More](#)

Parameter	Value
Product Name	CD10/MME Rabbit mAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43; CD10/MME
Species	Human
GeneID (Human)	4311
GeneID	4311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of human CD10/MME (NP_000893.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08473
Protein Weight	86kDa
Shipping	Ice bag