

Anti-Cd10/Mme Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08473

Item Number: CM0029566



Introduction

Rabbit polyclonal antibody against human CD10/MME (CALLA, Neprilysin) validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for cancer, neurodegeneration, and peptide hormone research.

[Learn More](#)

Parameter	Value
Product Name	CD10/MME Rabbit pAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43; CD10/MME
Species	Human
Gene ID (Human)	4311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-352 of human CD10/MME (NP_000893.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08473
Protein Weight	86kDa
Shipping	Ice bag