

Anti-Cd112 Rabbit Polyclonal Antibody For Wb, Elisa - P32507

Item Number: CM0021652



Introduction

Rabbit polyclonal antibody against mouse CD112 (Nectin-2). Validated for WB and ELISA, cross-reacts with human. Modulator of T-cell signaling and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	CD112 Rabbit pAb
Remarks/Alias	MPH; Pvr; Pvs; Cd112; Pvr12; nectin-2; CD112
Species	Mouse
Gene ID (Human)	5819
Gene ID (Mouse)	19294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-345 of mouse CD112 (NP_033016.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P32507
Protein Weight	57kDa
Shipping	Ice bag