


Anti-Cd115/Csf-1R Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P07333

Item Number: CM0029127



Introduction

High-quality rabbit polyclonal antibody targeting human CD115/CSF-1R (M-CSF receptor). Suitable for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with human, mouse, and rat. An essential tool for studies in macrophage biology, innate immunity, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD115/CSF-1R Rabbit pAb
Remarks/Alias	FMS; CSFR; FIM2; GP5C; HDLS; C-FMS; CD115; HDLS1; CSF-1R; BANDDOS; M-CSF-R; CD115/CSF-1R
Species	Human
GeneID (Human)	1436
GeneID	1436
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 873-972 of human CD115/CSF-1R (NP_005202.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07333
Protein Weight	108kDa
Shipping	Ice bag