

Anti-Cd116/Gm-Csfr α Rabbit Monoclonal Antibody For Wb, Elisa - P15509

Item Number: CM0006368



Introduction

Rabbit monoclonal antibody against human CD116/GM-CSFR α (CSF2RA, CDw116). Validated for WB and ELISA. Targets the low-affinity granulocyte-macrophage colony-stimulating factor receptor subunit alpha.

[Learn More](#)

Parameter	Value
Product Name	CD116/GM-CSFR α Rabbit mAb
Remarks/Alias	GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; alphaGMR; GMR-alpha; GMCSFR-alpha; GM-CSF-R-alpha; CD116/GM-CSFR α
Species	Human
GeneID (Human)	1438
GeneID	1438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-320 of human CD116/GM-CSFR α (NP_006131.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P15509
Protein Weight	27kDa/33kDa/38kDa/44kDa/46kDa/50kDa
Shipping	Ice bag