

Anti-Cd116/Gm-Csfra Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P15509

Item Number: CM0008225



Introduction

Recombinant rabbit monoclonal antibody against human CD116/GM-CSFR α , the GM-CSF receptor alpha subunit. Validated for Western Blot, Flow Cytometry, and ELISA applications with human reactivity. Multiple molecular weight isoforms detected: 27-50 kDa.

[Learn More](#)

Parameter	Value
Product Name	CD116/GM-CSFR α Rabbit mAb
Remarks/Alias	GMR; CD116; CSF2R; SMDP4; CDw116; CSF2RX; CSF2RY; GMCSFR; CSF2RAX; CSF2RAY; alphaGMR; GMR-alpha; GMCSFR-alpha; GM-CSF-R-alpha; CD116/GM-CSFR α
Species	Human
Gene ID (Human)	1438
Gene ID	1438
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-320 of human CD116/GM-CSFR α (NP_006131.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P15509
Protein Weight	27kDa/33kDa/38kDa/44kDa/46kDa/50kDa
Shipping	Ice bag