

Anti-Cd117 Rabbit Polyclonal Antibody For Wb, Elisa - P05532

Item Number: CM0024414



Introduction

Rabbit polyclonal antibody against mouse CD117/c-Kit (P05532), validated for WB and ELISA. Suitable for studying hematopoiesis, stem cell maintenance, and melanogenesis in mouse models.

[Learn More](#)

Parameter	Value
Product Name	CD117 Rabbit pAb
Remarks/Alias	W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3
Species	Human
GeneID (Human)	16590
GeneID	16590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117 (NP_001116205.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	109kDa
Shipping	Ice bag