

Anti-Cd117/C-Kit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10721

Item Number: CM0012380



Introduction

Rabbit monoclonal antibody against human CD117/c-Kit (c-Kit, CD117, SCFR) for Western blot, IF/ICC, IHC-P, and ELISA. Validated for human research in stem cell biology, oncology, and mast cell disorders.

[Learn More](#)

Parameter	Value
Product Name	CD117/c-Kit Rabbit mAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC
Species	Human
GeneID (Human)	3815
GeneID	3815
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 871-976 of human CD117/c-Kit (NP_000213.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag