

Anti-Cd117/C-Kit Polyclonal Antibody For Wb And Elisa - P10721

Item Number: CM0026713



Introduction

Rabbit polyclonal antibody against human CD117/c-Kit (KIT), validated for Western blot and ELISA. Reacts with human, mouse, rat. Suitable for stem cell, cancer, and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	CD117/c-Kit Rabbit pAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC; CD117/c-Kit
Species	Human
GeneID (Human)	3815
GeneID	3815
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 867-976 of human CD117/c-Kit (NP_000213.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag