

Anti-Cd117/C-Kit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P10721

Item Number: CM0014120



Introduction

Anti-CD117/c-Kit rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, rat. Recombinant protein immunogen. Research use in stem cell biology, hematopoiesis, melanogenesis, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	CD117/c-Kit Rabbit pAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC; CD117/c-Kit
Species	Human
GeneID (Human)	3815
GeneID	3815
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 210-412 of human CD117/c-Kit (NP_000213.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag