

Anti-Cd117/C-Kit Rabbit Monoclonal Antibody For Wb, Elisa - P05532

Item Number: CM0007073



Introduction

Rabbit monoclonal antibody targeting mouse CD117/c-Kit (stem cell factor receptor), validated for Western blot and ELISA. Suitable for studies in hematopoiesis, stem cell biology, and oncogenic signaling.

[Learn More](#)

Parameter	Value
Product Name	CD117/c-Kit Rabbit mAb
Remarks/Alias	W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3; CD117/c-Kit
Species	Mouse
GeneID (Human)	16590
GeneID	16590
Immunogen	Recombinant Protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	23KDa/109KDa
Shipping	Ice bag