

Anti-Cd117/C-Kit Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P05532

Item Number: CM0006661



Introduction

CD117/c-Kit Rabbit mAb validated for WB, FC, ELISA. Reacts specifically with mouse. Key tool for hematopoiesis, stem cell, and melanoma studies. SWISS: P05532.

[Learn More](#)

Parameter	Value
Product Name	CD117/c-Kit Rabbit mAb
Remarks/Alias	W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3; CD117/c-Kit
Species	Mouse
GeneID (Human)	3815
GeneID	16590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117/c-Kit (NP_001116205.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	109KDa/22KDa
Shipping	Ice bag