

Anti-Cd119 Polyclonal Antibody For Wb, Elisa - P15261

Item Number: CM0021176

Introduction

Rabbit polyclonal antibody targeting mouse CD119 (Interferon gamma receptor 1).
Validated for Western blot and ELISA, detects a 52 kDa membrane receptor essential for antimicrobial and antitumor immunity.

[Learn More](#)

Parameter	Specification
Product Name	CD119 Rabbit pAb
Remarks/Alias	Ifgr; CD119; Ifngr; Nktar; IFN-gammaR
Species	Mouse
Gene ID	15979
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-250 of mouse CD119. (NP_034641.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P15261
Protein Weight	52kDa
Shipping	Ice bag