


Anti-Cd11A Monoclonal Antibody For Wb, Fc - P20701

Item Number: CM0006320



Introduction

Anti-CD11a Rabbit Monoclonal Antibody validated for WB and FC in human. Recognizes Integrin alpha-L (LFA-1), a key immune adhesion receptor. Ideal for immunology and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	CD11a Rabbit mAb
Remarks/Alias	CD11A; LFA-1; LFA1A; CD11a
Species	Human
GeneID (Human)	3683
GeneID	3683
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a (NP_002200.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC
Cross Reactivity	Human
SWISS	P20701
Protein Weight	118kDa/129kDa/134kDa
Shipping	Ice bag