

Anti-Cd11A Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P24063

Item Number: CM0008430

Introduction

Mouse CD11a rabbit monoclonal antibody validated for Western blot, flow cytometry, and ELISA; cross-reacts with mouse and rat. The target is the 128 kDa integrin alpha-L chain (LFA-1) involved in leukocyte adhesion and immune functions.

[Learn More](#)

Parameter	Value
Product Name	CD11a Rabbit mAb
Remarks/Alias	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Species	Mouse
GeneID (Human)	728642
GeneID	16408
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 149-340 of mouse CD11a (NP_001240801.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P24063
Protein Weight	128kDa
Shipping	Ice bag