

Anti-Cd11A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P20701

Item Number: CM0007090



Introduction

Anti-CD11a rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human CD11a (ITGAL) with cross-reactivity to mouse and rat. Ideal for studying leukocyte adhesion and immune cell interactions.

[Learn More](#)

Parameter	Value
Product Name	CD11a Rabbit mAb
Remarks/Alias	ITGAL; CD11A; LFA-1; LFA1A; integrin alpha-L; CD11a
Species	Human
Gene ID (Human)	3683
Gene ID	3683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a (NP_002200.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20701
Protein Weight	119kDa/128kDa/134kDa
Shipping	Ice bag