


Anti-Cd11A/Lfa-1A/Itgal Monoclonal Antibody For Wb, Ihc-P, Elisa - P20701

Item Number: CM0002092



Introduction

Recombinant rabbit monoclonal antibody against human CD11a/LFA-1A/ITGAL (P20701), validated for WB, IHC-P, ELISA. Ideal for leukocyte adhesion, immune synapse, and phagocytosis research.

[Learn More](#)

Parameter	Value
Product Name	CD11a/LFA-1A/ITGAL Rabbit mAb
Remarks/Alias	CD11A; LFA-1; LFA1A; CD11a/LFA-1A/ITGAL
Species	Human
GeneID (Human)	3683
GeneID	3683
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1071-1170 of human CD11a/LFA-1A/ITGAL (P20701).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P20701
Protein Weight	129kDa
Shipping	Ice bag