

Anti-Cd11A/Lfa-1A/Itgal Rabbit Polyclonal Antibody For Wb, Elisa - P20701

Item Number: CM0021905



Introduction

Rabbit polyclonal antibody targeting human CD11a/LFA-1A/ITGAL, validated for WB and ELISA. Cross-reacts with human and rat. Suitable for studying leukocyte adhesion, immune cell interactions, and inflammatory signaling.

[Learn More](#)

Parameter	Value
Product Name	CD11a/LFA-1A/ITGAL Rabbit pAb
Remarks/Alias	CD11A; LFA-1; LFA1A; CD11a/LFA-1A/ITGAL
Species	Human
Gene ID (Human)	3683
Gene ID	3683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-306 of human CD11a/LFA-1A/CD11a/LFA-1A/ITGAL (NP_001107852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P20701
Protein Weight	129kDa
Shipping	Ice bag