

Anti-Cd11A/Lfa-1A/Itgal Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P20701

Item Number: CM0019341



Introduction

High-quality rabbit polyclonal antibody targeting human CD11a/LFA-1A/ITGAL. Validated for WB, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suited for leukocyte adhesion and immune cell studies.

[Learn More](#)

Parameter	Value
Product Name	CD11a/LFA-1A/ITGAL Rabbit pAb
Remarks/Alias	CD11A; LFA-1; LFA1A; CD11a/LFA-1A/ITGAL
Species	Human
GeneID (Human)	3683
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-340 of human CD11a/LFA-1A/ITGAL (NP_001107852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20701
Protein Weight	129kDa
Shipping	Ice bag