

Anti-Cd11B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11215

Item Number: CM0006617



Introduction

CD11b Rabbit mAb validated for Western Blot (WB), IHC-P, and ELISA. Suitable for human samples, the antibody targets integrin alpha-M (ITGAM) involved in leukocyte adhesion and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	CD11b Rabbit mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
GeneID (Human)	3684
GeneID	3684
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD11b (NP_000623.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag