

Anti-Cd11C Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P20702

Item Number: CM0011092



Introduction

CD11c Rabbit PolymAb® is a rabbit monoclonal antibody against human integrin alpha-X (CD11c), validated for WB, IF-P, IHC-P, ELISA. It reacts with human and cynomolgus monkey, and is ideal for studying monocyte adhesion and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	CD11c Rabbit PolymAb®
Remarks/Alias	CD11C; SLEB6
Species	Human
GeneID (Human)	3687
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 311-410 of human CD11c (NP_000878.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Cynomolgus monkey
SWISS	P20702
Protein Weight	128 kDa
Shipping	Ice bag