

Anti-Cd123 Rabbit Monoclonal Antibody For Fc, Elisa - P26952

Item Number: CM0004193

Introduction

CD123 Rabbit Monoclonal Antibody (IL-3R α) for flow cytometry and ELISA applications. Reacts with mouse CD123; recombinant immunogen corresponding to extracellular domain. Ideal for hematopoietic and cytokine signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD123 Rabbit mAb
Remarks/Alias	CD123; SUT-1; CDw123
Species	Mouse
GeneID (Human)	3563
GeneID (Mouse)	16188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-331 of mouse CD123 (NP_032395.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Mouse
SWISS	P26952
Protein Weight	43kDa
Shipping	Ice bag