

Anti-Cd127/IL-7R α Monoclonal Antibody () For Flow Cytometry - P16871

Item Number: CM0004717



Introduction

Recombinant rabbit monoclonal antibody against human CD127/IL-7R α , validated for flow cytometry. Recognizes human and cynomolgus monkey IL-7R α . Ideal for immunology, immune cell development, and disease target studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human/Monkey CD127/IL-7R α mAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha
Species	Human
GeneID (Human)	3575
GeneID	3575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-236 human CD127/IL-7R α (NP_002176.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag