

Anti-Cd127/IL7R Monoclonal Antibody For Wb, If-P, Elisa - P16871

Item Number: CM0003543



Introduction

Recombinant Anti-CD127/IL7R rabbit monoclonal antibody validated for WB, IF-P, ELISA; reactive with human, mouse, rat CD127. Ideal for T cell development and cytokine signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD127/IL7R Rabbit mAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha; CD127/IL7R
Species	Human
GeneID (Human)	3575
GeneID	3575
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human CD127/IL7R (P16871).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag