

Anti-Cd127/IL7R Rabbit Polyclonal Antibody For Wb, Elisa - P16871

Item Number: CM0016320



Introduction

High-quality rabbit polyclonal antibody against human CD127 (IL-7 receptor alpha), validated for Western blot and ELISA, cross-reacts with mouse. Useful for immunology and receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD127/IL7R Rabbit pAb
Remarks/Alias	ILRA; CD127; IL7RA; CDW127; IMD104; sIL-7R; Inc-IL7R; IL7Ralpha; IL-7Ralpha; IL-7R-alpha; CD127/IL7R
Species	Human
GeneID	3575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of human CD127/IL7R (NP_002176.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16871
Protein Weight	52kDa
Shipping	Ice bag