

Anti-Cd130/Gp130 Rabbit Monoclonal Antibody For Flow Cytometry - Q00560

Item Number: CM0005084



Introduction

High-specificity rabbit monoclonal antibody against mouse CD130/gp130, validated for flow cytometry. Target molecular weight 102kDa. Suitable for IL-6 signaling pathway research and cytokine receptor studies.

[Learn More](#)

Parameter	Value
Product Name	CD130/gp130 Rabbit mAb
Remarks/Alias	CD130; gp130; D13Ert699e; 5133400A03Rik; CD130/gp130
Species	Mouse
Gene ID (Human)	3572
Gene ID (Mouse)	16195
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00560
Protein Weight	102kDa
Shipping	Ice bag