


Anti-Cd132 Rabbit Polyclonal Antibody For Wb, Elisa - P34902

Item Number: CM0021940



Introduction

Rabbit polyclonal antibody against mouse CD132 (IL-2RG, gammaC), validated for WB and ELISA. Recognizes rat CD132. Suitable for cytokine signaling and immune research studies.

[Learn More](#)

Parameter	Value
Product Name	CD132 Rabbit pAb
Remarks/Alias	gc; p64; [g]c; CD132; gamma(c)
Species	Mouse
Gene ID	16186
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-263 of mouse CD132 (NP_038591.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P34902
Protein Weight	42kDa
Shipping	Ice bag