


Anti-Cd132/IL-2 R Gamma Polyclonal Antibody For Wb And Elisa - P31785

Item Number: CM0023829



Introduction

Rabbit polyclonal antibody targeting human CD132/IL-2 R gamma (common gamma chain). Validated for Western blot and ELISA; cross-reacts with mouse. Suitable for studying cytokine receptor signaling and primary immunodeficiency research.

[Learn More](#)

Parameter	Specification
Product Name	CD132/IL-2 R gamma Rabbit pAb
Remarks/Alias	P64; CIDX; IMD4; CD132; SCIDX; IL-2RG; SCIDX1; CD132/IL-2 R gamma
Species	Human
GeneID (Human)	3561
GeneID	3561
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-369 of human CD132/IL-2 R gamma (NP_000197.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31785
Protein Weight	42kDa
Shipping	Ice bag