

Anti-Cd133 (Prom1) Rabbit Monoclonal Antibody For Wb, Ip, If-P, Elisa - O54990

Item Number: CM0005871

Introduction

CD133 rabbit monoclonal antibody validated for WB, IP, IF-P, ELISA. Reacts with mouse and rat. Detects ~97 kDa Prominin-1, a key regulator of cell differentiation and retinal disk morphogenesis.

[Learn More](#)

Parameter	Value
Product Name	CD133 Rabbit mAb
Remarks/Alias	Prom; AC133; CD133; Prom-1; Proml1; 4932416E19Rik
Species	Mouse
GeneID (human)	8842
GeneID	19126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-107 of mouse CD133 (NP_001157049).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O54990
Protein Weight	97 kDa
Shipping	Ice bag