

Anti-Cd134 Rabbit Polyclonal Antibody For Wb, Elisa - P47741

Item Number: CM0019703

Introduction

Rabbit polyclonal antibody targeting mouse CD134/OX40, validated for WB and ELISA, with cross-reactivity to rat. Ideal for T-cell costimulatory pathway and long-term immunity studies.

[Learn More](#)

Parameter	Value
Product Name	CD134 Rabbit pAb
Remarks/Alias	Ox40; ACT35; CD134; Ly-70; Txgp1; TXGP1L
Species	Mouse
GeneID	22163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-211 of mouse CD134 (NP_035789.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P47741
Protein Weight	30kDa
Shipping	Ice bag