


Anti-Cd137 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P20334

Item Number: CM0007744



Introduction

Mouse CD137 (4-1BB) Rabbit Monoclonal Antibody validated for WB, IF/ICC, and ELISA. Recognizes the T-cell costimulatory receptor TNFRSF9 to support immuno-oncology and T-cell biology research.

[Learn More](#)

Parameter	Value
Product Name	CD137 Rabbit mAb
Remarks/Alias	ILA; Ly63; 4-1BB; Cd137; CDw137; A930040I11Rik
Species	Mouse
GeneID (Human)	3604
GeneID	21942
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	P20334
Protein Weight	28kDa
Shipping	Ice bag