

Anti-Cd138 Monoclonal Antibody For Flow Cytometry - K7Gnf5

Item Number: CM0005491

Introduction

High-quality -conjugated Rabbit anti-Pig CD138 monoclonal antibody (SDC1, syndecan-1) validated for flow cytometry in pig. Recombinant protein immunogen. Ideal for studying exosome biogenesis, cell-matrix interactions, and dental development. 32kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Pig CD138 mAb
Remarks/Alias	SDC1, syndecan-1
Species	Pig
GeneID (Human)	6382
GeneID	100519770
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Pig
SWISS	K7GNF5
Protein Weight	32kDa
Shipping	Ice bag