

Anti-Cd138/Syndecan-1 Rabbit Monoclonal Antibody For Wb - P18827

Item Number: CM0009802



Introduction

High-quality rabbit monoclonal antibody targeting human CD138/Syndecan-1 (SWISS: P18827) validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for IVD research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	CD138/Syndecan-1 Rabbit mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan; CD138/Syndecan-1
Species	Human
Gene ID (Human)	6382
Gene ID	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-254 of human Syndecan 1 (P18827)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag