

Anti-Cd13/Anpep Monoclonal Antibody For Fc - P15144

Item Number: CM0004751



Introduction

Rabbit monoclonal antibody targeting human CD13/ANPEP, specifically designed for flow cytometry (FC) applications. Immunogen spans amino acids 69-967. Useful for immunology and virology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD13/ANPEP mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
Species	Human
GeneID (Human)	290
GeneID	290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag