

Anti-Cd13/Anpep Monoclonal Antibody For Flow Cytometry - P15144

Item Number: CM0007994



Introduction

rabbit monoclonal antibody against human CD13/ANPEP (P15144), validated for flow cytometry. Targets a broad-specificity aminopeptidase involved in peptide digestion, angiogenesis, and viral entry. Ideal for immunophenotyping and functional studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD13/ANPEP mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
Species	Human
GeneID (Human)	290
GeneID	290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag