

Anti-Cd13/Anpep Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P15144

Item Number: CM0003107



Introduction

High-quality monoclonal antibody targeting human CD13/ANPEP, validated for WB, IF-P, IHC-P, ELISA applications. Cross-reacts with mouse and rat. Ideal for studies on aminopeptidase N function and host-virus interactions.

[Learn More](#)

Parameter	Value
Product Name	CD13/ANPEP Rabbit mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150; CD13/ANPEP
Species	Human
Gene ID (Human)	290
Gene ID	290
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human CD13/ANPEP (P15144)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag