

Anti-Cd13/Anpep Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - P15144

Item Number: CM0003503



Introduction

CD13/ANPEP Rabbit monoclonal antibody validated for WB, IHC-P, FC, ELISA; immunogen recombinant human CD13 (aa 69-967); 110kDa; for diagnostic and research use targeting aminopeptidase N.

[Learn More](#)

Parameter	Value
Product Name	CD13/ANPEP Rabbit mAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150; CD13/ANPEP
Species	Human
GeneID (human)	290
GeneID	290
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag