

Anti-Cd13/Anpep Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P15144

Item Number: CM0030363



Introduction

CD13/ANPEP Rabbit pAb is a polyclonal antibody targeting aminopeptidase N. Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Useful for studying protein digestion, angiogenesis, and viral infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	CD13/ANPEP Rabbit pAb
Remarks/Alias	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150; CD13/ANPEP
Species	Human
GeneID (Human)	290
GeneID	290
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 668-967 of human CD13/ANPEP (NP_001141.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15144
Protein Weight	110kDa
Shipping	Ice bag