


Anti-Cd14 Monoclonal Antibody For If/Icc, Ihc, Elisa - P08571

Item Number: CM0003506



Introduction

Rabbit monoclonal antibody against human CD14, suitable for immunofluorescence, immunohistochemistry, and ELISA. A key tool for studying innate immunity, LPS signaling, and monocyte/macrophage biology.

[Learn More](#)

Parameter	Value
Product Name	CD14 Rabbit mAb
Remarks/Alias	CD14
Species	Human
GeneID	929
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human CD14 (P08571)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag