


Anti-Cd14 Rabbit Monoclonal Antibody For Wb - P10810

Item Number: CM0009848



Introduction

Mouse CD14 rabbit monoclonal antibody validated for WB, FC, ELISA. Reacts with mouse samples, detects ~39 kDa protein. Suitable for innate immunity and LPS signaling studies.

[Learn More](#)

Parameter	Value
Product Name	CD14 Rabbit mAb
Remarks/Alias	CD14
Species	Mouse
Gene ID (Human)	929
Gene ID (Mouse)	12475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-345 of mouse CD14 (NP_033971.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse
SWISS	P10810
Protein Weight	39k Da
Shipping	Ice bag