


Anti-Cd141/Thrombomodulin Polyclonal Antibody For Wb, Ihc-P, If/Icc, And Elisa - P07204

Item Number: CM0029856



Introduction

Rabbit polyclonal thrombomodulin antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat CD141/THBD; regulates coagulation and inflammation. Ideal for cardiovascular and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	CD141/Thrombomodulin Rabbit pAb
Remarks/Alias	TM, THRM, AHUS6, BDCA3, CD141, BDCA-3, THPH12, CD141/Thrombomodulin
Species	Human
GeneID (Human)	7056
GeneID	7056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-515 of human CD141/Thrombomodulin (NP_000352.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07204
Protein Weight	60kDa
Shipping	Ice bag