


Anti-Cd141/Thrombomodulin Rabbit Monoclonal Antibody For IHC-P, IF/ICC, Fc, Elisa - P07204

Item Number: CM0002353



Introduction

Rabbit monoclonal antibody targeting human CD141/Thrombomodulin (THBD), validated for IHC-P, IF/ICC, flow cytometry, and ELISA. Suitable for coagulation, inflammation, and endothelial cell biology research.

[Learn More](#)

Parameter	Value
Product Name	CD141/Thrombomodulin Rabbit mAb
Remarks/Alias	TM; THRM; AHUS6; BDCA3; CD141; BDCA-3; THPH12; CD141/Thrombomodulin
Species	Human
GeneID (Human)	7056
GeneID	7056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-515 of human CD141/Thrombomodulin (NP_000352.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P07204
Protein Weight	60kDa
Shipping	Ice bag