


Anti-Cd141/Thrombomodulin Rabbit Polymab® For Wb, Ihc-P, If/Icc, Fc, Elisa - P07204

Item Number: CM0006516



Introduction

High-specificity rabbit monoclonal antibody against human CD141/Thrombomodulin. Validated for WB, IHC-P, IF/ICC, FC, ELISA. Recombinant immunogen. Ideal for coagulation, inflammation, and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	CD141/Thrombomodulin Rabbit PolymAb®
Remarks/Alias	TM; THRM; AHUS6; BDCA3; CD141; BDCA-3; THPH12
Species	Human
GenelD (human)	7056
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-515 of human CD141/Thrombomodulin. (NP_000352.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P07204
Protein Weight	60kDa
Shipping	Ice bag