

## Anti-Cd142 (Tissue Factor) Monoclonal Antibody For If/Icc, Fc - P13726

Item Number: CM0011610



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD142/Tissue Factor, validated for immunofluorescence and flow cytometry. Recombinant immunogen, cross-reacts with human. Ideal for hemostasis and coagulation research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD142 mAb                                                                                        |
| Remarks/Alias    | TF; TFA; CD142                                                                                                     |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 2152                                                                                                               |
| Gene ID          | 2152                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 34-251 of human CD142 (NP_001984.1). |
| Source           | Rabbit                                                                                                             |
| Category         | Monoclonal Antibodies                                                                                              |
| Application      | IF/ICC, FC                                                                                                         |
| Cross Reactivity | Human                                                                                                              |
| SWISS            | P13726                                                                                                             |
| Protein Weight   | 33kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |